

Compact
Powerful
Precise



ColorFlex[®] L2 Tomato

*Level UP your Color Confidence with the new
ColorFlex L2 Tomato Spectrophotometer.*

As simple as you want, as powerful as you need.

Elevate your color accuracy and user experience like never before. Compact and ergonomic build requires minimal space while offering enhanced durability and reliability. Designed to streamline and expedite your quality control processes, ColorFlex L2 Tomato Spectrophotometer is expertly crafted to measure the color of tomatoes in all their processed forms, including paste/purée, sauce, ketchup, juice, and fresh tomatoes.




Explore the Next-Gen ColorFlex®...L2 Tomato



#1 Color Measurement Solution for Tomatoes

HunterLab Tomato Reference Standard: HunterLab, the creators of the original tomato color reference standard, offers the only reference tile that traces back to the discontinued BCR and USDA standards, ensuring the most accurate red color measurements.




All the Tomato Scores and Indices You Need

Measure Tomato in Every Form: Tomato Paste Score (TPS), Tomato Sauce Score (TSS), Tomato Catsup Score (TCS), Tomato Juice Score (TJS), a/b ratio, Fresh Tomato Color Index (FTCI) & Lycopene index.

45°a:0° Annular Illuminated Optical Geometry

Superior Consistency and Color Accuracy: Enhancing the efficiency of Quality Control workflows and delivery of superior quality products repeatedly.



Compact with Sealed Spillproof Case

Space-Saving and Easy Maintenance: Occupies minimal space and is designed to withstand the demands of busy lab environments.



Tomato Paste
Tomatoes

Tomato Sauce
TOMATO KETCHUP

tomato juice
Tomato Soup




Powerful Standalone Functionality

No computer required: With huge onboard processing and storage capacity, ColorFlex L2 Tomato can function as an independent color workstation. Connect monitor, keyboard, and mouse if desired.




Modern Tablet User Interface with Touchscreen Color Display

Durable and Efficient: Facilitates accurate color measurement in seconds, with a robust industrial glass design suitable for harsh environments.




New & Improved Color Software

Essentials 2.0: Take precise and reliable color measurement in seconds with next-gen user-friendliness and seamlessly analyze multiple tristimulus and spectral data sets. It's that powerful.



Self-Guided and Intuitive Out-of-Box Startup

Time-Saving and User-Friendly: Ready to use straight out of the box with an integrated wizard for quick training, reducing operator time and easing concerns related to staff training and turnover.



Built-in Camera for Sample Viewing and Screen Recording

Precision and Communication: Ensures accurate sample positioning to eliminate measurement discrepancies and allows easy sharing of recordings for swift quality control decisions.



Tomato Paste
Tomatoes


Tomato Sauce
TOMATO KETCHUP

tomato juice
Tomato Soup



Built-in Advanced Diagnostics

On-Board Diagnostic Health Check: Ensures the HunterLab ColorFlex L2 Tomato is always operating at peak performance, continuing to build on the legendary history of reliability of the ColorFlex family of instruments.



Advanced Data Export and Communication

Full suite of data communication options: Features USB ports, HDMI, and support for mouse and keyboard, enabling easy data export/printing.



Connectivity with LIMS and SPC systems

Simplifies QC workflows: Aids in trend analysis and workflow simplification across manufacturing sites.

CONTACT US

We've worked with thousands of companies, but you'll feel like our one and only. If color accuracy is crucial to your business, you've come to the right place—with the right people. Let's get in touch. Ask a question, request a quote or just say hello!

OUR OFFICES & DISTRIBUTION CENTERS

HunterLab
 11491 Sunset Hills Road
 Reston, VA 20190-5280
 Website: www.hunterlab.com
 Email: sales@hunterlab.com
 Ph: 703-471-6870
 Fax: 703-471-4237

Corporate Headquarters

HunterLab Europe
 Christian Jansen, Managing Director
 Dr.-August-Einsele-Ring 15
 82418 Murnau, Germany
 Website: www.hunterlab.com
 Email: Christian.Jansen@hunterlab.com
 Ph: +49 8841 9464

European Office

HunterLab Asia-Pacific
 Heng Su Kiang, APAC Sales Manager
 28-2, Plaza KLTS
 No. 99, Jalan Gombak
 Setapak
 53000 Wilayah Persekutuan Kuala Lumpur
 Malaysia
 Website: www.hunterlab.com
 Email: sales-asiapacific@hunterlab.com
 Ph: +(60) 340 324 749

Asia-Pacific Offices

HunterLab Asia-Pacific
 Leo Cheung, Regional Manager
 Flat E, 13F Leahander Centre
 28 Wang Wo Tsai Street Tsuen Wan
 New Territories Hong Kong
 Website: www.hunterlab.com
 Email: sales-asiapacific@hunterlab.com
 Ph: +852-2406 1723

For more information, please contact HunterLab at 703-471-6870, sales@hunterlab.com or visit www.hunterlab.com